TANITA

Dual Frequency Body Composition Monitor

Model: RD-445

en

Instruction manual

Read this Instruction Manual carefully and keep for future reference.

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en**Table of contents**

Introduction	1
Features and Functions	1
For Your Safety	2
Preparations Before Use	4
Automatic Shut-Down Function	4
Setting the Date and Time	5
Initial Setup	6
Taking Measurements	7
How To See Your Segmental Measurements	9
Recall Function	9
Taking Measurements Without Registering (Guest User Measurements / Weight Only Measurements)	10
Troubleshooting	10
Specifications	11

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Introduction

Thank you for purchasing this TANITA Body Composition Monitor.

This monitor is one of a wide range of home health care products produced by TANITA.

This Instruction Manual will guide you through the setup procedures and outline the key features of this monitor. Additional information on healthy living can be found on our website www.tanita.eu. TANITA products incorporate the latest clinical research and technological innovations.

All data is monitored by the TANITA Medical Advisory Board to ensure accuracy.

Safety precautions

Do not use the Body Composition Monitor if you have an electronic medical implant such as a pacemaker, as it passes a low-level electrical signal through the body which may interfere with the operation of the implant.

Pregnant women should only use the weight function. Other functions are not intended for use when pregnant.

Important notes for users

This Body Composition Monitor is intended for adults and children (age 5-17) who have inactive to moderately active lifestyles, and adults with athletic lifestyles.

Thanks to major investment in the latest BIA Technology and sports science research, TANITA has upgraded the Athlete Mode function to make measurements more accurate and suitable for a wider range of users.

Who should use Athlete Mode?

Adults aged 18 years and over who either:

- Train or exercise for 12 hours or more a week and have been doing so for at least six months.
- Are body builders.
- Are professional athletes who want to monitor their progress at home.
- Have a lifetime history of fitness and who used to do more than 12 hours a week but do less now.

Recorded data may be lost if the unit is used incorrectly or is exposed to electrical power surges.



TANITA cannot accept any responsibility for losses incurred due to the loss of recorded data.

TANITA cannot accept any responsibility for damages or losses incurred through the use of this product or any third-party claims.


Note: Body fat percentage estimates will vary with the amount of water in the body, and can be affected by dehydration or over hydration due to such factors as alcohol consumption, menstruation, illness, intense exercise, etc.

Features and functions

Measuring Platform

AA-size batteries

For your safety

This section explains precautionary measures to be taken to avoid injury to the users of this device and others, and to prevent damage to property. Please familiarize yourself with this information to ensure safe operation of this equipment.

Warning

Failure to follow instructions highlighted with this mark could result in death or severe injury.

Caution

Failure to follow instructions highlighted with this mark could result in injury or damage to property.

This mark indicates actions that are prohibited.

This mark indicates instructions that must always be followed.

Warning

Prohibited

This monitor must not be used by people with pacemakers or other medical implants.

This monitor passes a weak electrical current through the body which could interfere with and cause the malfunction of electrical medical implants, resulting in serious harm to the user.

Interpretation of measurements (e.g. evaluation of measurements and formulation of exercise programmes based on the measurements) must be performed by a professional.

Implementing weight loss measures and exercise programmes based on self-analysis could be harmful to your health. Always follow the advice of a qualified professional.

Do not step on the edge of the platform. Do not jump onto, or on the monitor.

Doing so may cause you to fall or slip, resulting in serious injury.

Do not place this monitor on a slippery surface such as a wet floor.

Doing so may cause you to fall or slip, resulting in serious injury. It may also cause water to get inside the monitor resulting in damage to internal components.

Caution

Prohibited

Do not stand on the display screen or control panel.

Do not insert your fingers into gaps or holes.

Doing so may result in injury.

Do not use rechargeable batteries. Do not use old batteries together with new batteries, or a mix of different types of batteries at the same time.

Doing so may cause battery fluid to leak or the batteries to become excessively hot and rupture, resulting in damage to the monitor or injury.

Never allow an infant or small child to carry the monitor.

Store out of reach of small children.

The child may drop the monitor, resulting in injury.

Small children may become caught up in the cables, resulting in injury.

This monitor is intended for home use only. This monitor is not intended for professional use including hospitals, medical or fitness facilities. It is not equipped with the components required for heavy usage. Using the monitor in a professional location will invalidate the warranty.

Do not pull out the cable unless using the monitor for measurement, as doing so may damage the monitor.

Do not pull out the cable beyond the red mark. Do not pull the cable too forcefully.

If the cable does not automatically return into the monitor unit, untwist the cable, pull it out approximately 20cm and then try returning it again.

Place the monitor on a hard, flat surface where there is minimal vibration to ensure safe and accurate measurement.

Avoid excessive impact or vibration to the monitor.

Doing so may damage the monitor or cause it to malfunction, or may cause setting to be erased.

Do not use or store this monitor in the following locations.

- Locations that are exposed to direct sunlight, near to heaters, or exposed directly to air from air conditioning units, etc.
- Locations where there are severe temperature changes or a large amount of vibration.
- Locations with high levels of humidity or moisture.

Doing so may damage the monitor or cause it to malfunction.

Do not step on the monitor when your body is wet (immediately after a shower or bath, etc.)

Be sure to thoroughly dry your body and feet before using the monitor. Failure to do so may cause you to slip or fall, resulting in injury. Also, the monitor cannot take accurate measurements if your body or feet are wet.

Keep away from water.

Do not use boiling water, benzene, thinner, etc. to clean the monitor.

If the monitor becomes dirty, soak a soft cloth in water or neutral household detergent, wring it out thoroughly and wipe the monitor clean, then use a dry cloth to wipe the monitor dry.

Do not attempt to disassemble the monitor. Do not modify this equipment in any way.

Do not wash the monitor in water.

Do not submerge this monitor in water.

Doing so may damage the monitor or cause it to malfunction.

Be sure to insert the batteries with the polarities (+, -) in the correct position.

If the batteries are incorrectly inserted, battery fluid may leak and damage the floor. (If you do not intend to use this monitor for a long period of time, it is advisable to remove the batteries before storage.)

Do not use the monitor near other products such as microwave ovens that emit electromagnetic waves.

Doing so may cause the monitor to malfunction, preventing accurate measurements.

Another person should assist persons with disabilities who may not be able to take a measurement alone.

Avoid measuring after over-eating or over-drinking, or when severely dehydrated.

This may cause inaccurate measurements. For greater accuracy, avoid using directly after waking up. Use at the same time of day each time, at least two hours after the last meal.

Avoid taking measurements after spending time in a sauna, bath or after intense exercise.

This may cause inaccurate measurements. Please take measurements after sufficient rest.

Avoid taking measurements when your physical condition is poor (hangover, diarrhea, fever, etc.).

Doing so may cause error in the measured value. Resume measuring once you have recovered.

Perform measurement with as little clothing on as possible.

Weight of clothing worn will be included in measurements.

Be sure that your toes and heels are correctly aligned with the electrodes on the measurement platform.

The body fat percentage displayed may be low or result in an "Error" (error) display if contact between the electrodes and the soles of the feet is poor.

If you do not undress, always remove your socks or stockings, and be sure the soles of your feet are clean before stepping on the measuring platform.

Use in a stable location. Errors in measurements may occur if the device is used in an unstable location.

An accurate reading will not be possible if your knees are bent or you are in a sitting position.

Do not move during measurement.


This causes inaccurate measurements.

Ensure arms are not touching sides and inner thighs are not touching each other during measurement. If necessary, place a dry towel between arm and side and/or between thighs.


Always...

Preparations before use

Open as shown.

Direction of the batteries.

Inserting the Batteries

Open the battery cover on the back of the measuring platform. Insert the supplied AA-batteries as indicated.

Note:

- Be sure that the polarity of the batteries is set properly. If the batteries are incorrectly positioned the fluid may leak and damage floors.
- If you do not intend to use this unit for a long period of time, it is advisable to remove the batteries before storage.

Positioning the Scale

Place the measuring platform on a hard, flat surface where there is minimal vibration to ensure safe and accurate measurement.


Note: To avoid possible injury, do not step on the edge of the platform.

Do not pull the Hand-electrodes cable beyond the red band.

Changing the Weighing Unit

Turn on the scale, then press the button on the back of the scale to switch the weighing unit between kg, lb, and st-lb.


Set Your Location Code

Different locations across the world have different gravitational pulls. These differences affect the accuracy of the weight measurement. You can obtain the most accurate weight readings by selecting the correct gravity setting on the monitor according to your geographical location.

1. When the power is off, press and hold the USER button.
2. Click [here](#) or scan the QR code to view your location code, and press the Up or Down button to select your location code.
3. Press the SET button to enter and save the setting.

***Note:** Only at the time of purchase, the location setting screen will appear without having to press the USER button.

Automatic shut-down function

The automatic shut-down function shuts off the power automatically in the following cases:


- If you interrupt the measurement process. The power will shut down automatically within 20 to 60 seconds.
- During programming, if you do not touch any of the buttons within 60 seconds.*
- After you have completed the measuring process within 60 seconds.

***Note:** If the power shuts off automatically, repeat the steps from "Initial Setup" (page 6).

Setting the date and time


1.Set the year

Insert batteries into the scale and then press the ON/OFF button to start the Year setting. Press the Up or Down button to set year and then press the SET button. The unit will beep once to confirm.


2.Set the month

Press the Up or Down button to set month and then press the SET button. The unit will beep once to confirm.


3.Set the day

Press the Up or Down button to set day and then press the SET button. The unit will beep once to confirm.

4.Set the hour


Press the Up or Down button to set hour and then press the SET button. The unit will beep once to confirm.

5.Set the minute

Press the Up or Down button to set minute and then press the SET button.

The unit will beep twice to confirm the programming. The power will then shut down automatically after 20-60 seconds.

Note:

- After you replace the batteries or if you want to reset date and time, press and hold the Down button to begin the date and time programming process. Replacing the batteries will cause the loss of date and time settings.
- 60 seconds or more of inactivity will cause automatic shut off during this process. No data will be saved as a result.
- If you make a mistake or wish to turn the unit off before you complete the programming, press the ON/OFF key to turn the unit off. No data will be saved as a result.
- Only after turning on the device for the first time; the date and time setting screen will appear after the location setting.

Initial setup

The unit can be operated only if data has been programmed into one of the personal data memories.

1. Press the ON/OFF key while the Grip (Hand-electrodes) is mounted on the platform scale

"0.0" will be shown.

2. Start setting

Press the SET button.

The unit will beep to confirm activation, the Personal data number will be displayed, and the display will flash.

Note:

- If you don't operate the unit for 60 seconds after the unit has been turned on, the unit will turn off automatically.
- If you make a mistake or want to turn the unit off before you have finished programming it, press the ON/OFF key to force quit.

3. Select a Personal data Number

Press the Up or Down button to select a Personal data number (1, 2, 3, 4 and 5).

Once you reach the personal data number you wish to use, press the SET button.

The unit will beep once to confirm.

4. Set your Birthday

Set the date of your birthday (Year / Month / Day) using the Up or Down button.

Press the SET button to confirm at each programming.

The unit will beep once to confirm.

Note: The range of age is from 5 – 99.

5. Select Female or Male

Use the Up or Down button to scroll through Female (♀), Male (♂), Female/Athlete (♀♂) and Male/Athlete (♂♂) settings, then press the SET button.

The unit will beep once to confirm.

6. Specify the Height

The display defaults to 170 cm (5' 7.0") (range of height is from 90 cm – 220 cm or 3' 0.0" – 7' 3.0"). Use the Up or Down button to specify Height and then press the SET button. The unit will beep once to confirm. The display shows the personal data one time to confirm, and then the display shows "0.0" and is ready for measurement.

7. Measuring your weight and body fat

Hold the Hand-electrodes and step on to the scale after "0.0" is shown on display within 30 seconds.

Do not step off until the unit beeps five times.

The readings will be shown for 40 seconds.

The unit will then shut off and the readings will be stored in the memory.












Note: If you step onto the scale before "0.0" appears in the display, then "Error" message will be displayed, and you have to step off and on again. If you do not step onto the scale within about 30 seconds after "0.0" being displayed, the scale will be shut off automatically. Please try the setting and storing personal data process again.

Note:

If you change the birthday and/or gender, the unit will show "CLr".

If you select "YES" using the Up or Down button, the existing personal data will be deleted.

If you select "no" using the Up or Down button, no new data will be saved (existing data will be unchanged) then the unit will shut off.

Taking Measurements

After programming your personal data you are ready to take a reading.

1. Press the ON/OFF button while the Grip (Hand-electrodes) is mounted on the platform scale

"0.0" will be shown.


2. Get your readings

Step onto the scale and then pick up the Grip (Hand-electrodes).

Your Body weight will be shown first.


Continue to stand on the scale.

Your personal profile is shown on the display and then the Body weight, the BMI and the Body fat% reading are shown for 40 seconds.

Note: Do not step off until the unit beeps twice.

Make sure all of your fingers contact electrodes.

Ensure that your arms are full extended and your elbows do not touch your body.

Make sure your feet contact all of the electrodes.

Do not bend your knees.

Do not pull the hand-electrodes cable beyond the red band.


Note:

If the wrong personal number appears on the display, you can change the personal number by using the USER button.

Press the USER button until your personal number and profile appears.



After you select your correct personal number and press the SET button, the display will show your correct readings automatically.

The readings will be stored in your personal number memory.

Continued on the next page.

Immediately after measuring, you can press the Result button while still standing on the scale to view your results. You can also step off the scale and then press the Result button to view your results again.


Note: For children (age 5-17), the unit will only display Body weight, BMI, Body fat%, Bone mass and Muscle mass.

How to see your segmental measurements?

What does segmental measurements mean?


A segmental measurement provides a detailed analysis of your muscle mass and fat percentage in your arms, legs, and torso. Using both foot and hand electrodes, the scale measures each arm, leg, and torso separately for a more accurate assessment.

Press the Muscle Mass or Body Fat% button repeatedly to view the values for the right arm, left arm, right leg, left leg, and torso.

Recall function

- Press the ON/OFF key. "0.0" will be shown.
- Press the USER button until your personal number appears.
- Press the SET button and PAST button to see the previous readings.

Press the  button or the  button to show previous readings of each segment/per body part when the previous readings are displayed on the screen.

Note:

For children (age 5-17), the unit will only display the Body weight, the BMI and the Body fat%.


Taking Measurements Without Registering

(Guest User Measurements/
Weight-Only Measurements)

Guest User Measurements

Guest mode allows you to programme the unit for a one-time use without resetting a Personal data number.

- Press the ON/OFF button while the Hand-electrodes is mounted on the platform scale.
- "0.0" will be shown.
- Hold the Hand-electrodes and the USER button until "Guest" appears.
- Press the SET button and input your personal information such as Age, Sex, Height. Please see page 6 for more details.
- After "0.0" will be shown, step on the scale holding the Hand-electrodes.


Do not step off until the unit beeps twice.

The readings will be shown for 40 seconds.

Weight-Only Measurements

- Press the ON/OFF button while the Grip (Hand-electrodes) is mounted on the platform scale.
- After 2 or 3 seconds, "0.0" will appear in the display and the unit will beep.
- When "0.0" appears in the display, step onto the platform. The Body weight value will appear.
- The power will then shut down automatically.

Troubleshooting

If the following problems occur ...then...

- **A wrong weight format appears followed by kg, lb, or st-lb**

Refer to "SWITCHING THE WEIGHT MODE" (page 4).

- ** appears on the display, or all the data appears and immediately disappears.**

Batteries are low. When this message appears, be sure to replace the batteries immediately, since weak batteries will affect the accuracy of your measurements. Change all the batteries at the same time with new batteries.

Note:

Your settings and past memory will not be erased from the memory when you remove the batteries.

- ** appears while measuring.**

Are the personal data settings correct?

Did you select the wrong personal data number?

Please stand on the platform keeping movement to a minimum.

The unit cannot accurately measure your weight if it detects movement.

Measurement items have exceeded the measurement range. (Measurements cannot be taken)

- **The Body Fat Percentage reading does not appear.**

Make sure socks or stockings are removed, and the soles of your feet are clean and properly aligned with the guides on the measuring platform.

- ** appears while measuring.**

Reading cannot be obtained if the weight capacity is exceeded.

Specifications

RD-445

Weight Capacity

200kg/440lb/31st6lb

Weight Increments

0-100kg/0.05kg 100kg-200kg/0.1kg

0-200lb/0.1lb 200-440lb/0.2lb

0-15st0lb/0.1lb 15st0lb-31st6lb/0.2lb

Applicable Age Range

Body Mass Index

5-99 years old

Total & Segmental Body Fat (%)

5-99 years old

Body Water %

18-99 years old

Total & Segmental Muscle Mass

5-99 years old

Muscle Quality

18-99 years old

Physique Rating

18-99 years old

Bone Mass

5-99 years old

Basal Metabolic Rate

18-99 years old

Metabolic Age

18-99 years old

Visceral Fat Rating

18-99 years old

Personal Data

5 memories

Time accuracy

Within ± 1 minute per month (+18 - +28°C)

Power Supply

4×AA batteries

Power Consumption

60mA Maximum

Measuring Current

50kHz, 5kHz, 100μA

Disposal

This product is an electronic device. Please dispose of this product appropriately, not as general household waste. Be sure to follow the regulations in your area when disposing of this product.

Not allowed to mix batteries with consumer wastes!

As consumer you are legally bound to return used or discharged batteries. You can deposit your old batteries at the public collecting points in your town, or wherever the corresponding batteries are sold and specifically marked collecting boxes have been set up. In case of scrapping the apparatus, the batteries should be removed from it and deposited at the collecting points as well.

The device features radio interference suppression in compliance with EC regulation 2004/108/EC

Notice:

These marks can be found out on contaminant-containing batteries:

Pb Pb = Battery contains lead,

Cd Cd = Battery contains cadmium,

Hg Hg = Battery contains mercury.

- The product design and specifications may be changed at any time without prior notice.

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Federal Communications Commission and Canadian ICES Notice

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules and Canadian ICES-003. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio or television technician for help.

Modifications

The FCC requires the user to be notified that any changes or modifications made to this device that are not expressly approved by Tanita Corporation may void the user's authority to operate the equipment.

Representative TANITA Europe B.V.

Hoogoorddreef 56-E
1101 BE Amsterdam
The Netherlands
Tel: +31-20-560-2970
Fax: +31-20-560-2988
www.tanita.eu

TANITA Health Equipment H.K. Ltd.

Unit 301-303 3/F Wing On Plaza,
62 Mody Road, Tsimshatsui East,
Kowloon, Hong Kong
Tel: +852-2838-7111
Fax: +852-2838-8667
www.tanita.com

TANITA India Private Limited

Unit 502 of 'A' Wing,
Mittal Commercial,
Marol, Andheri East,
Mumbai 400059 INDIA
TEL: +91 22 3192 6107
info@tanita.co.in

Manufacturer TANITA Corporation

1-14-2, Maeno-Cho, Itabashi-ku,
Tokyo, Japan 174-8630
Tel: +81-3-3968-7048
Fax: +81-3-3968-2661
www.tanita.co.jp
ISO 9001 Certified

TANITA Corporation of America, Inc.

2625 South Clearbrook Drive
Arlington Heights, IL 60005 USA.
Tel: +1-847-640-9241
Fax: +1-847-640-9261
www.tanita.com